


Tri*nt*rac®

INTRAORAL DISTRACTORS



CE 1984

Triintra[®]

INTRAORAL DISTRACTORS

TITANIUM FIXATION PLATE TYPES



FP- H- 060 – L/R Titanium Fixation plate "H" Type, 6 holes



FP- ST- 040 Titanium Fixation plate, Straight Type, 4 holes



FP- LF1– L/R Titanium Fixation plate, Lefort 1 Type, 6 holes



FP- ST- 030- L/R Titanium Fixation plate, Lefort 1 Type, Straight, 3 holes



FP- T- 060 – L/R Titanium Fixation plate, "T" Type, 6 holes



FP- DS- 080 Titanium Fixation plate, Double Side Straight Type, 8 holes

Tri[®]intrac

INTRAORAL DISTRACTORS

DISTRACTION CYLINDERS



DC- 0020 Distraction Cylinder Ø 4.0mm X 20 mm



DC- 0015 Distraction Cylinder Ø 4.0mm X 15 mm



DC- 0010 Distraction Cylinder Ø 4.0mm X 10 mm

FLEXIBLE DRIVER



FW- 070 Flexible driver 70mm



FW- 040 Flexible driver 40mm



FW- 030 Flexible driver 30mm

T-WRENCH



TW-010 Distraction Wrench

61.SCR.5.0 Titanium Micro Screw Ø 1.6mm X 5.0mm

61.SCR.6.0 Titanium Micro Screw Ø 1.6mm X 6.0mm

61.SCR.7.0 Titanium Micro Screw Ø 1.6mm X 7.0mm

61.SCR.9.0 Titanium Micro Screw Ø 1.6mm X 9.0mm



TWIST BONE DRILL BIT



61.DR.60/25 Stainless Steel Bone Drill Bit ø 1.2mm X 25mm X 60mm

Triintrac®

INTRAORAL DISTRACTORS

SCREW DRIVER




D7.750 Stainless Steel Screw Driver

TROCAR SET



E1.180 Trocar Set



MANUFACTURER: UCMEDE TIBBİ CİHAZLAR LTD. ŞTİ.
GERSAN SANAYİ SİTESİ 2307. Sok. No: 22 Ergazi - ANKARA - TURKEY
TEL: +90 312 256 28 38 - 255 73 26 FAX: +90 312 255 73 25
web: www.trimed.com.tr e-mail: trimed@trimed.com.tr